Instruction for use BC-San 100

BC-San 100

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Inhalt:
12 x 100 ml e

Contents

1.	User instructions	4
2.	Symbols used on product information, labels and instructions for use	4
3.	Target group	5
4.	Intended purpose – intended use	5
5.	Virucidal final disinfection of surfaces of medical devices	6
5.1.	Pretreatment	6
5.2.	Cleaning and pre-disinfection	6
5.3.	Preparation for virucidal final disinfection	6
5.4.	Virucidal final disinfection	6
6.	Processing of wipe dispenser boxes	7
6.1.	Preparation	7
6.2.	Cleaning and disinfection	7
6.3.	Drying	7
6.4.	Filling	7
6.5.	Information label „LOT and EXPIRY“	7
7.	Processing of procedural water supply bottles/tanks	8-9
7.1.	Cleaning and disinfection of the bottle holding fixture	8-9
7.2.	Cleaning and disinfection of procedural water supply bottles	8-9
7.3.	Contact time	8-9
8.	Active ingredient basis	10
9.	Efficacy	10
10.	Order data	10
11.	Hazard and precautionary statements	11
12.	Disposal instructions	11
13.	Report of incidents	11

1. User instructions

Please attend to the following instructions in order to guarantee a trouble-free, economic and safe work.

2. Symbols used on product information, labels and instructions for use

	temperature limitation (here: storage and transportation temperature)
	Use before (date of expiry) YYYY-MM
	Batch number XXXXXX
	Please follow the instruction for use
	EEC-symbol for finished packs
	CE-symbol with Notified Body (0123) for Medical Devices of classes IIa, IIb or III
	Order number
	Manufacturer, Company ALPRO MEDICAL GMBH, St. Georgen
	Keep product away from sunlight

3. Target group

This instruction for use is intended for all users (doctors, dentists and professionals in charge of the processing of medical devices as well as hygiene professionals) in medical and dental facilities.

—> **Application only by medical professionals**

4. Intended purpose – intended use

BC-San 100 CE⁰¹²³

Solution for the virucidal final disinfection of surfaces of medical devices, such as light guides, polymerisation lamps, heads of intraoral cameras/scanners, as well as for the disinfection/processing of wipe dispenser boxes and pressure bottles (1 - 1.5 L) of stand-alone procedural water supply systems such as the ALPRO-BCS (BottleCareSystem).

5. Virucidal final disinfection of surfaces of medical devices

5.1. Pretreatment

Remove visible coarse contamination by means of a cellulose wipe.


5.2. Cleaning and pre-disinfection

Wipe the surface of the medical devices by means of a TB-effective disinfection wipe and let it take effect according to the specified contact time; e.g. with PlastiSept eco Wipes or MinutenWipes.

5.3. Preparation for virucidal final disinfection

Disconnect (if possible) the instrument's part to be disinfected with BC-San 100 from the "base".

Soak a suitable cellulose wipe with BC-San 100.

5.4. Virucidal final disinfection

Wipe thoroughly by means of the BC-San 100 soaked cellulose wipe the instrument sectors to be disinfected and let it take effect for at least 5 minutes.

6. Processing of wipe dispenser boxes

6.1. Preparation

Open the wipe dispenser box and dispose the remaining disinfection wipes and disinfection solution after the expiry.

6.2. Cleaning and disinfection

Wipe thoroughly all internal parts of the open wipe dispenser box, incl. the tear-off element and the closing cap, with a suitable cellulose wipe soaked with BC-San 100.

6.3. Drying

Leave recipient and lid open to dry for 10 minutes.

6.4. Filling

Put new refill bag in the processed wipe dispenser box and cut open along the marked line.

→ **Note:** Please refer for further steps to the instruction for use of e.g. MaxiWipes reel.

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ALPRO
ALPRO MEDICAL GMBH

7. Processing of procedural water supply bottles/tanks

7.1. Cleaning and disinfection of the bottle holding fixture

Disconnect the procedural water supply bottle from the bottle holding fixture. Generally, clean and disinfect per working day at end of the day (e. g. by means of MinutenWipes or with suitable cellulose wipe soaked with undiluted BC-San 100) the bottle holding fixture and the riser tube (see pictures).

7.2. Cleaning and disinfection of procedural water supply bottles

For the treatment of the internal surfaces (e. g. of BCS-bottles) pour the content of BC-San 100 in the empty procedural water supply bottles/tanks (see pictures). Afterwards, fill the procedural water supply bottles with drinking water up to the screwing (screw thread/quick fastening).

7.3. Contact time

Let the solution in the procedural water supply bottles take effect at least 6 hours (e. g. overnight).

Keep the filled procedural water supply bottles during disinfection closed in the preparation room or mount it as usual at the holding fixture of the procedural water supply system.

Attention: Do not put pressure on the procedural water supply bottles during the disinfection on the ALPRO-BCS (BottleCareSystem)!

- Dispose the solution with BC-San 100 in the sink after the contact time resp. before reconnection of the water supply system.
- Afterwards, rinse the procedural water supply bottle with water.

8. Active ingredient basis:

Sodium hypochlorite solution with 1-2% active chlorine.

9. Efficacy:

- Sporicidal (prEN 17126, Bacillus subtilis) in 10 minutes
- Virucidal (EN 16777, Parvo viruses) in 5 minutes

Test criteria: clean conditions, room temperature

10. Order data BC-San 100

REF	Description	
3199	12 x 100 ml bottle in cardboard box, incl. instruction for use	

11. Hazard and precautionary statements

Warning. Causes skin irritation. Causes serious eye irritation. Harmful to aquatic life with long lasting effects.

Wear protective gloves/eye protection/face protection.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

IF EYE IRRITATION PERSISTS: Get medical advice/attention.

Not suitable for unalloyed metals.

Contact with acids liberates toxic gas.

Protect from heat and direct sunlight.

12. Disposal instructions

The ready-to-use solution can be disposed of in highly diluted form via the waste water. The undiluted product or product residues must be disposed of as hazardous waste in compliance with the Directive 2008/98/EC on waste as well as national and regional regulations. Do not dispose of the undiluted product via the waste water. Leave product in the original container as possible. Do not mix with other waste materials. Packaging contaminated with product is considered as hazardous waste and must be disposed of accordingly.

13. Report of incidents

Please report any serious incident that may occur in connection with the product to:

ALPRO MEDICAL GMBH

E-Mail: quality@alpro-medical.com

Tel.: +49 7725 9392-0

and the national authorities

Information on the preparation with regard to the work steps and all current work instructions, safety data sheets, instructions for use etc. can be found on our website:

www.alpro-medical.com

If you have any questions, please do not hesitate to contact us at **www.alpro-medical.com** or by phone **+49 7725 9392-0**. You are also welcome to contact our sales representatives by phone or e-mail.