



Efficacy

- ▶ bactericidal incl. MRSA and TBC
- ▶ yeasticidal
- ▶ virucidal against enveloped viruses in Carrier Test DVV 2012 incl. HBV, HCV, HIV, herpes simplex, influenza (incl. H1N1, H5N1), BVDV, vaccinia, corona viruses (e. g. MERS-CoV) as well as non-enveloped noro viruses

Application

- ▶ completely moisten surfaces and objects to be disinfected
- ▶ let the disinfectant take effect according to the spectrum of activity
- ▶ wipe off with a disposable wipe (e. g. PlastiSept-Wipes eco)

Available as

- | | |
|----------|--|
| REF 4645 | 500 ml foam bottle |
| REF 4646 | 5 l canister |
| REF 4647 | 5 l cubitainer |
| REF 4630 | box with 120 wipes |
| REF 4631 | carton with 20 x 120
wipes refill |
| REF 4632 | box with 200 jumbo
wipes |
| REF 4633 | carton with 12 x 200
jumbo wipes refill |



ALPRO MEDICAL GMBH

ALPRO MEDICAL GMBH • Mooswiesenstr. 9
78112 St. Georgen • Germany
Tel: 0 77 25 / 93 92 0 • Fax: 0 77 25 / 93 92 91
www.alpro-medical.de • info@alpro-medical.de

REF 80350250GB

The NOVELTY of ALPRO



PlastiSept eco

**Free of QAC,
alcohol and
aldehydes,
TBC effective**

ALPRO
ALPRO MEDICAL GMBH

Gentle surface disinfection

PlastiSept eco

PlastiSept eco is a ready for use disinfectant especially suitable for alcohol sensitive surfaces of medical devices and instruments.

Advantages

- ▶ without substances damaging the material
- ▶ without allergenic or sensitizing substances
- ▶ particularly gentle on all materials
- ▶ free of alcohol, aldehydes and QAC

Field of application

- ▶ surfaces of medical devices, e. g. treatment units, upholstery, head rests and suspension tables
- ▶ surfaces of instruments, e. g. turbines, contra-angles and ultrasonic scalers
- ▶ surfaces in the X-ray area, e. g. chin rests at OPG
- ▶ surfaces of acrylic glass (also for solarium beds)

