



Do you have rust on your instruments?

Are you tired of manual cleaning of your endo files and surgical burs?

ALPROZYME alkaline, enzymatic granulates can help you to solve these problems.

Before AlproZyme



After AlproZyme



Contact time

15-30 min in soaking bath

5-15 min in ultrasonic bath

Max contact time recommended 60 min.

1 sachet makes 2L of the solution!

Designed for soaking baths, ultrasonic baths and thermal washers, AlproZyme prevents corrosion and dissolves the inorganic and organic debris. Your instruments will look brand new!

Reorder directly from Henry Schein Halas (reference # RD-3753M) and receive 10% OFF on your first order.