


Adapters for WL-series

as at: 8 July 2015



REF 3521 **Adapter 01** (plastic),
Universal adapter with Luer connector for transmission instruments (ISO 9168) and media channels of the medical driven motors. It's used for the application and the blowing out of cleaning solutions (WL-clean, WL-cid, maintenance oil). Also suitable for blowing out by using the multifunction syringe.
Note: With this adapter all medical hollow part instruments equipped with a "Luer connector" can be precleaned, e. g. with WL-clean.

Adapters for spray water and spray air channels*:



REF 3522 **Adapter 02** (stainless steel),
for spray water and spray air channels
Adapter for spray water and spray air channels of all **hand pieces and contra-angles with ISO-interface**.



REF 3523 **Adapter 03** (stainless steel),
for spray water and spray air channels
Adapter for spray water and spray air channels for **turbines with "Sirona"-interface**.



REF 3524 **Adapter 04** (stainless steel),
for spray water and spray air channels
Adapter for spray water and spray air channels for **turbines with „Multiflex“- (KaVo/SciCan) interface**.



REF 3526 **Adapter 06** (stainless steel),
for spray water and spray air channels
Adapter for spray water and spray air channels for **turbines with „Roto-Quick“-interface (W&H)**.




REF 3530 **Adapter 07** (plastic),
for spray water and spray air channels
Adapter for long shank contra-angles solely suitable for spray water and spray air channels of **Sirona T1 Classic-Line**.

*note the current approvals of instrument manufacturers!

All adapters indicated above are not suitable for maintenance with oil!

Adapters for complete internal cleaning*:




REF 3532 **Adapter 02/G** (stainless steel)
for complete internal cleaning (gearing and media channels)
Adapter for **contra-angles and hand pieces with ISO-interface**
like adapter 02 but for the complete internal cleaning/disinfection*
(gearing incl. head of instrument).



REF 3533 **Adapter 03/G** (stainless steel)
for complete internal cleaning (gearing and media channels)
Adapter for **turbines with "Sirona"-interface** like adapter 03
but for the complete internal cleaning/disinfection*
(gearing incl. head of turbine).



REF 3534 **Adapter 04/G** (stainless steel)
for complete internal cleaning (gearing and media channels)
Adapter for **turbines with „Multiflex“-(KaVo/SciCan) interface**
like adapter 04 but for the complete internal cleaning/disinfection*
(gearing incl. head of turbine).



REF 3535 **Adapter 05/G** (stainless steel)
for complete internal cleaning (gearing and media channels)
Adapter for **turbines with „Bien-Air“-interface** for the complete
internal cleaning/disinfection*(gearing incl. head of turbine).



REF 3536 **Adapter 06/G** (stainless steel)
for complete internal cleaning (gearing and media channels)
Adapter for **turbines with „Roto-Quick“- interface (W&H)**
like adapter 06 but for the complete internal cleaning* with
WL-clean/ WL-Blow (gearing incl. head of turbine).



REF 3539 **Adapter 09/G** (stainless steel)
for complete internal cleaning (gearing and media channels)
Adapter for **turbines with "NSK"-interface** for the complete
internal cleaning/disinfection* (gearing incl. head of turbine).

*note the current approvals of instrument manufacturers!

All adapters indicated above are not suitable for maintenance with oil!

Adapters for WL-series

as at: 8 July 2015



Morita

REF 3561 **Adapter 10/G** (stainless steel)
for complete internal cleaning (gearing and media channels)
Adapter for **turbines with "Morita"-interface** for the
complete internal cleaning/disinfection*
(gearing incl. head of turbine).



Acteon

REF 3563 **Adapter 11/G** (stainless steel)
for complete internal cleaning (gearing and media channels)
Adapter for **ACTEON-AIR N GO Airpolisher**
"Supra/Perio nozzle" for the complete internal cleaning*
with WL-clean, WL-cid, WL-Blow.



REF 3571 **Adapter U-1 + U-1a Kombi** (stainless steel)
Combined adapter
Combined adapter to clean* surgical contra-angles and hand
pieces externally cooled with the WL-series. Consisting of U-1
rinsing adapter for ultrasonic tips and U-1a screw attachment
for hose plug connection.

*note the current approvals of instrument manufacturers!



REF 3541 **Introductory set**
for complete internal cleaning
(Adapter 02/G, 04/G) incl. 1 replacement set of o-rings
Introductory set with 2 adapters

REF 3528	O-Ring-Set for adapter 02 3 x 2 seals	REF 3529	O-Ring-Set for adapter 03 und 06 3 x 3 seals
REF 3543	O-Ring-Set for adapter 02/G 1 x 2 seals	REF 3545	O-Ring-Set for adapter 03/G 1 x 2 seals
REF 3544	O-Ring-Set for adapter 04/G 1 x 3 seals	REF 3546	O-Ring-Set for adapter 06/G 1 x 1 seal

All adapters indicated above are not suitable for maintenance with oil!

Adapters for ultrasonic and airscaler tips*:



REF 3570 **Adapter U-1** (stainless steel)
rinsing adapter for ultrasonic tips M3 „internal thread“
from EMS, Mectron, KaVo, Komet



REF 3572 **Adapter U-2** (stainless steel)
rinsing adapter for ultrasonic tips M3 „external thread“
from EMS, Mectron, KaVo, Komet, NSK




REF 3573 **Adapter U-3** (stainless steel)
rinsing adapter for ultrasonic tips M3/a „internal thread“
from Satelec, NSK



REF 3574 **Adapter U-4-S** (stainless steel)
rinsing adapter for ultrasonic tips in the Piezosurgery „internal thread“
from NSK (VarioSurg), Mectron (Piezosurgery).



REF 3575 **Adapter U-5-S** (stainless steel)
rinsing adapter for ultrasonic tips in the ultrasonic surgery „external thread“ from Sirona (Siroprep and Siroimplant).



Note the manufacturers approvals!

All adapters indicated above are not suitable for maintenance with oil!